



## Instructions for Shipment from Australia

1. We recommend using our internet-based DHL platform, the Express Premium Portal (EPP-Access), to order DHL-Pickup and to create the corresponding Shipment Waybill: [www.arminlabs.com/dhl](http://www.arminlabs.com/dhl)  
Please **print out your individual DHL Waybill** at the end of the shipment order and place it behind the transparent window of the DHL plastic bag (step 4).
2. Alternatively, **we will gladly take care of the creation of the DHL Express Shipment Waybill**. In this case, please contact us early enough (at least 3-4 days in advance; mind the time difference) to guarantee a smooth shipment. We need your name, the pick-up address and your telephone number. Afterwards, you will receive a waybill created by ArminLabs via e-mail. Please call your local DHL Express Service to book the collection of your shipment. You can find your local DHL service number next to your country here: <http://www.arminlabs.com/content/forms/phone-support-dhl-arminlabs.en.pdf>  
Please tell the DHL service employee our account number **956 801 388**.

1		<ul style="list-style-type: none"> <li>• Provide the patient's name on every tube</li> <li>• Take the blood samples</li> <li>• Put the tubes into the Styrofoam box, close it with the rubber band and put it into the airtight zip lock bag</li> <li>• Complete and sign the <b>Test Order Form</b></li> <li>• Put the Styrofoam box and the Test Order Form into the carton and close it</li> </ul>
2		<ul style="list-style-type: none"> <li>• Fill in the green dashed boxes of the <b>Confirmation Letter TDI</b> ("your name", "waybill number below barcode", "day of shipment") and <b>3 Proforma Invoices</b> Sign the <b>3 Proforma Invoices</b></li> </ul>
3		<ul style="list-style-type: none"> <li>• Put the completed <b>Confirmation Letter TDI</b> and the <b>3 Proforma Invoices</b> behind the transparent window of the DHL plastic bag.</li> </ul>

4	 <p>The image shows a DHL shipping label for a medical specimen. It includes the sender's address (Musterfirma, Helsinki, Finland), the recipient's address (Arminlabs GmbH, Augsburg, Germany), and shipping details such as 'DE-MUC-MUC', 'RET', and '0.3 kg 1/1'. There are three barcodes at the bottom.</p>	<ul style="list-style-type: none"> <li>Place the <b>Shipment Waybill</b> (either the one you created in the DHL online platform or the one you received from ArminLabs) behind the transparent window of the DHL plastic bag in front of the Confirmation Letter TDI and the 3 Proforma Invoices.</li> </ul>
5	 <p>The image shows a small cardboard box with a 'SARSTEDT' logo being placed into a clear DHL plastic shipping bag.</p>	<ul style="list-style-type: none"> <li>Put the filled carton into the provided <b>transparent DHL plastic bag</b></li> </ul>
6	 <p>The image shows the DHL plastic bag from the previous step being sealed with a piece of adhesive tape.</p>	<ul style="list-style-type: none"> <li>Close the DHL plastic bag with the help of the adhesive tape on the bag</li> </ul>
7	<p>Your <b>Test-Kit</b> is ready for shipping with DHL</p>	

We will gladly help you if you have any questions:

Phone: 0049 821 78093150  
E-mail: [info@arminlabs.com](mailto:info@arminlabs.com)